

**To all
SWISS-TRAC
users**

CH-6403 Küssnacht, 20.04.2015/j.-

SWISS-TRAC Wheelchair towing machine

Dear SWISS-TRAC user,

we can confirm the following:

- All SWISS-TRAC are fitted with air- and seaworthy lead batteries that do not contain liquid acids. **All batteries comply with the ICAO special provision A67.** The SWISS-TRAC can be safely transported with any vehicle (airplane, train, ship or car). The battery must not be disassembled or removed from its housing.
- The wheelchair towing machine is a mobility aid that belongs to your wheelchair. It weighs approx. 65 kg and measures roughly 80x50x58 cm (LxWxH).
- All airlines are obliged to carry mobility aids for the disabled free of charge, including a **SWISS-TRAC.**
- If your airline applies an additional charge, we advise you to make your reservation with another large airline company instead.

We advise you to inform your airline **that the SWISS-TRAC must never be transported on the luggage conveyor belt.**

Advise: Please use for air transport our SWISS-TRAC transport banner.

We wish you a safe journey and a great holiday!

Yours Sincerely,

ATEC Engineering Ltd.



J. Jakober Ing. FH
Managing Director