

ATEC Ingenieurbüro AG
 Breitenstrasse 1
 6403 Küssnacht
 Switzerland

Tel. +41 41 854 80 20
 Fax +41 41 854 80 21
 E-Mail info@swisstrac.ch
 MWST-Nr. CHE-106.299.803

SWISS•TRAC®

www.swisstrac.ch

**To the attention of
 Swiss-Trac users
 and transport companies**

CH-6403 Küssnacht, 09.01.2026

Official Transport Certificate and Battery Confirmation from the Manufacturer of Swiss-Trac Wheelchair Add-on Drives

Dear Swiss-Trac user,
 dear transport company

The manufacturer confirms the following:

- All wheelchair tractor devices of the brand **Swiss-Trac** can be safely transported within any vehicle like **airplane, train, ship, public transport (e.g. bus, metro), car and others**.
 The batteries must not be disassembled nor removed. No further precaution is necessary.
- The Swiss-Trac is equipped with non-spillable AGM lead batteries (WCBD, IATA DGR 2.3.2.2).
They are air- and seaworthy and comply with the IATA/ICAO special provision A67.
The 300 Wh limit refers to lithium batteries, excluding AGM batteries.
- The wheelchair tractor is a mobility aid that belongs to the wheelchair.
 Weight mod. SWT-1/SWT-1S/SWT-1HD: approx. 70 kg / 155 lbs
 Weight mod. SWT-2: approx. 55 kg / 122 lbs resp. see technical specification.

Addition for airlines:

- Airlines are obliged to carry mobility aids for people with disabilities free of charge, incl. Swiss-Trac.
- If your airline applies an additional charge, we advise you to make your reservation with another airline company instead.
- Inform your airline that the **Swiss-Trac must never be transported on the luggage belt**.
- Please use our [Swiss-Trac transport banner](#) for air transport.

Addition for public transport companies:

- The Swiss-Trac is tilt-proof.
- The person in the wheelchair has more stability whilst the Swiss-Trac is connected to the wheelchair.

We wish you a safe and pleasant journey.

Yours sincerely,



ATEC Ing. Büro AG
 Breitenstrasse 1
CH-6403 Küssnacht
 Tel. +41 (0)41 854 80 20

J. Jakober, Ing. FH, Managing Director
ATEC Ingenieurbüro AG